



## PRODUCT INSTRUCTION SHEET

REFERENCE      Orbital Seater with Knee Support (805)

### PURPOSE

Reduces strain on carer and enables safe and comfortable transfers for people who require additional support when standing.

### DESCRIPTION

High tensile steel revolving transfer platform. Two revolving laser cut steel discs with a top frame bolted to the discs and support handles welded to the top frame. 2 x 100mm Dia rubber tyres allow manoeuvrability. Polyurethane knee support (Code 806PUSUPPORT) (350mm x 160mm) clamps onto rear back rails.

### TECHNICAL SPECIFICATIONS

|   |                  |
|---|------------------|
| Width between arms:                         | 17" (430mm)      |
| Height of front handles from platform base: | 30" (760mm)      |
| Height of top bar from platform base:       | 39" (990mm)      |
| Width:                                      | 26" (660mm)      |
| Depth:                                      | 27" (685mm)      |
| Disc diameter:                              | 18" (450mm)      |
| Product weight:                             | 17kg             |
| Recommended user weight:                    | 35 Stone (222kg) |

### USAGE REQUIREMENTS

1. Check castors are securely fitted before during and after use. This is important since the lock nuts holding the two castors in place also secure the axle /tilt bar.
2. The unit is designed to be used on LEVEL SURFACES only.
3. Do not attempt to move the unit from one position to another with a user actually standing on it.
4. Please ensure that the locking device is fully located before transferring a user from one sitting position to another.
5. Periodically check that the top frame is always fully secured to the revolving base plate.
6. Always ensure that the edges of the non-slip surface, on top of the turning disc are fully flush with the surface of the disc.
7. Frames require an occasional wipe over with a soft cloth-use general purpose neutral detergent and warm water. Do not use harsh abrasive materials. To disinfect, clean as above and wipe over with a solution of 1000 parts per million chlorine disinfectant (haztabs/precept) rinse and dry.

## PRODUCT INSTRUCTION SHEET

### **INTENDED PURPOSE:**

The Secure Turn Patient Turner has been designed to rotate a standing patient during a transfer from e.g. wheelchair to a chair.

### **PRECAUTIONS:**

Ensure these instructions are fully read and understood before use and please keep in a safe place for future reference.

Please note the patient must not be transported using the device.

Do not exceed the maximum user weight limit of 35st (222Kg).

Ensure the handle is facing the wheel and that the brake is locked before manoeuvring the Secure Turn.

Never assist a patient to alight from the wheeled side of the Secure Turn.

Transfers should only be carried out to the left or right through 45° or 90°.

Patient alighting should only be performed when the brake is engaged and the Secure Turn is stabilised by a carer.

The Secure Turn should be inspected on a regular basis for signs of damage or wear that might affect safe use.

Do not clean this device using an industrial cleaner or high-pressure hose as this may affect the smooth operation of the turntable.

### **INSTRUCTIONS:**

1. Familiarise yourself with the operation of the Secure Turn before attempting to transfer a patient, in particular the brake locking mechanism
2. To manoeuvre the Secure Turn into position, with the handle facing the wheels and brake locked in position, tilt the device back onto the wheels and move into position.
3. Place the seat/chair that the patient is transferring to into position at 45° or 90° to the patient's current position.
4. Place the turntable up to the patient's feet with the wheels facing the carer.
5. The patient should position their feet, evenly spaced, on the turntable mat. The patient may need to move to the edge of the seat before rising to reduce the effort required.
6. When the patient needs to raise themselves from a sitting position. e.g. from below the height of the handles of the secure turn, then the carer must stabilize the unit by placing their foot on the axle (midway between the 2 wheels) to prevent tilt. Alternatively, if the patient needs to lower themselves onto the secure turn e.g. from a bed, then the carer should stabilize the unit by placing one foot on either of the lateral stabilizing lugs to prevent tilt.
7. With the Carer stabilising the Secure Turn, the patient raises themselves into a standing position.
8. When the patient is standing and stable (holding on to one of the pairs of handles on the frame), the carer releases the brake and rotates the turntable to the required position allowing the brake to lock into position.
9. The carer stabilises the Secure Turn whilst the patient lowers themselves onto the seat/chair.
10. After the patient has removed their feet from the turntable, manoeuvre the Secure Turn away from the patient.

**PRODUCT INSTRUCTION SHEET**

**800 Orbital Transfer Platform – Instructions for use**



**Care & Cleaning**

Wipe clean with a soft cloth, moistened with disinfectant solution.

**WARRANTY**

Your Drive Devilbiss branded product is warranted to be free of defects in materials and workmanship for one year from date of purchase. This device was built to exacting standards and carefully inspected prior to shipment. This warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the customer. In the event of a defect covered by this warranty, we will at our option supply parts or replace the device. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, castors and grips, which are subject to normal wear and need periodic replacement. If you have any queries or require full warranty conditions, please contact your Drive stockist. This does not affect your statutory rights.