




# The Ultima

TRANSFER BENCH WITH SLIDING  
SEAT AND FOLD UP LEGS

---

PRODUCT CODE  
**BA06042**

# ASSEMBLY INSTRUCTIONS



1. Unfold the shower chair by opening the leg set toward both sides, the leg set is securely fixed when a “click” sound heard.
2. Fold the shower chair by releasing the levers inward (see below drawing marked direction) on the leg set till the bolt break away from the connecting holes.
3. Adjust the extension tubes of each leg set so that the bench seat is level and at a comfortable height for the user. Push the spring button in and slide the leg tube to the next position.
4. Snap in the bench back; align the plugs and holes and push downwards.
5. Set the hand hold into position on the tub side of the bench, align the plugs and the holes then tighten the screws from below the bench.
6. Soap caddy can be used on either side of the bench. Align the tab and push into the centre hole.



## OPERATION

1. Place transfer bench in tub so the suction cup feet are **INSIDE** the tub. (If the preference is for the seat to face the opposite way, legs can be switched to the other side of the bench by depressing brass push pins, removing legs and inserting them into the opposite leg receptacles.)
2. Adjust height of legs to desired height by depressing brass push pins and setting to proper height. Note: Legs should be set to ensure seat is level, taking into account that the legs inside the tub will have to be set to a lower height setting because of the difference in height from most tubs to floors.
3. To move seat, pull the red seat positioning lever located under the seat away from the frame and slide seat until it clicks into place. Note: there are two (2) positions where seat will be secured to frame, one on the left side and one on the right side.

## FOLDING

To fold for storage or transport, stand bench vertically on one set of legs. On one side of leg assembly, pull both red locking levers away from frame simultaneously until both locking pins disengage. Fold legs down. Turn transfer bench over and repeat for second set of legs. Remove backrest, handle and soap dish if necessary.

## SPECIFICATIONS



1	Overall Height	815-915mm
2	Seat Depth	470mm
3	Suction Cup Leg width	450mm
4	Overall Width	1030mm
5	Rubber Feet Leg Width	440mm
6	Seat Height	485-590mm



## **CAUTION**

When entering or exiting shower or bath via the transfer bench, make sure the seat is in one of the locked positions.

## **WEIGHT CAPACITY**

300 LBS (136 KG)

## **WARRANTY**

Your CareCo branded product is warranted to be free of defects in materials and workmanship for one year from the date of purchase.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your CareCo device or this warranty, please contact your dealer.

