

CareCo

Comfort Sleep Washable Bed Pads

Product Code: BA05047



Maximum Coverage
Breathable,
Hypoallergenic &
Antimicrobial fabric




Non-Slip Backing
Absorbs up to 2 litres
per 1 metre²



Machine Washable
Machine washable
for more than
300 cycles



Locks in Liquids
100% water and
leak proof



Soft Quilted: 80% cotton
20% polyester fabric

High Density: 80%
polyester, 20% rayon
non-woven fabric

Super Absorbent:
Breathable & waterproof
polyurethane membrane

Anti-Slip: 100% polyester
knitted fabric

- Extra large size for total coverage and greater absorption
- Mattress and sheets stay dry and stain-free all night
- Lay it on top of your bottom sheet

34in x 36in
86cm x 91cm

Distributed by CareCo
Made in China

