

# Eden

High Back Chair

PRODUCT CODE

RR05081

MAXIMUM WEIGHT

**152kg**

DIMENSIONS

**75 cm** WIDE

**80 cm** DEEP

**117 cm** HIGH

SEAT WIDTH


**49 cm**

SEAT DEPTH

**45 cm**

SEAT HEIGHT

**39 - 52 cm**



- Unique air cell in seat for enhanced pressure care
- Features soft and comfortable 4-way stretch material
- 4-way stretch is breathable, water resistant and fire retardant
- Covers can be easily replaced
- Includes rear leg castors for transporting chair
- Telescopic leg height adjustment, seat depth adjustment
- Ergonomic back support

**WARRANTY:** Products are to be fit for purpose and of excellent quality and performance. For valid warranty claims CareCo will, at their discretion, replace / repair / refund items mutually agreed to be defective.

This chair is covered by a 1 year manufacturers warranty that warrants the goods will be free from material defects from manufacturing

This warranty does not cover items such as arm pads and seat cushions that are subject to wear and tear

Any damage or defect of any nature occurring from the misuse, abuse of the product, improper operation or improper storage is not covered. The warranty starts from the date of receipt of your product.