

# MOBILIS

FURNITURE RASIERERS




The 3" & 5" Mobilis Furniture Raisers provide you with a simple solution to get in and out of beds and chairs with ease and comfort.

The smooth edged design of the Furniture Rasiers mean they minimise scuffs and dents to your furniture and hardware. Made of heavy-duty PP material, the raisers are built for maximum durability.

- Suitable for legs that are between 2cm and 8cm in diameter.
- Maximum weight - 600kg (when using all 4 raisers)

## 3" Mobilis Furniture Raiser: DL06268.3




Height: **97mm**

Bottom Width: **130mm**

Inside depth: **17mm**

Inside width: **73mm**

## 5" Mobilis Furniture Raiser: DL06268.5



Height: **152mm**

Bottom Width: **164mm**

Inside depth: **21mm**

Inside width: **69mm**