



## **TriGlide** Cruiser



UP TO  
**35 MILE RANGE**



**AUTOMATIC SAFE-BRAKING SYSTEM**



**FULLY ADJUSTABLE  
COCKPIT**



**COMFORT & STYLE**



**TWIN DIAL  
DIGITAL DASH**



***The TriGlide Cruiser Mobility Scooter is the perfect blend of performance and comfort.***

With a top speed of 8mph, you're ready for the roads. The throttle lever lets you dial down the speed for travelling on pavements when you need to and the automatic safe-braking system keeps you in control.



## TriGlide Cruiser

### SPECIFICATIONS

Height	115cm
Width	76cm
Length	166cm
Max Speed	8mph
Wheel Size	14"
Full Suspension	Yes
Product Weight	195kg
Battery Size	85ah   100ah (12v x 2)
Max User Weight	28.6st
Control Panel Cover	Yes
Digital Display	Yes
Pannier Box Storage	Yes

